Pulsossimetro MANUALE D'USO A3/MAN-002 V2.4

www.healforce.com

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Shenzhen Hexin Zondan Medical Equipment Co. Ltd.
Floor 14,Block D,Dianlian Technology Building,the Crossing between South Circle
Road and South Full road,Guangming District
518106 Shenzhen
PEOPLE'S REPUBLIC OF CHINA

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9 051884 250021

Declaration of Conformity

Manufacturer Shenzhen Hexin ZONDAN Medical Equipment Co., Ltd.
Address Floor 14,Block D,Dianlian Technology Building,the Crossing between South Circle Road and South Fuli Road ,Guangming District,Shenzhen,China
European Representative Shanghai International Holding Corp. GmbH (Europe)
Product Eiffestraße 80 20537 Hamburg GERMANY
Model Code Pulse Oximeter
UMDNS A2\A3\A4\A5
17148

Classification

The medical device has been classified to class II a according to rule 10 of Annex IX of the Directive 93/42/EEC.

Conformity Assessment Route: Annex II without chapter 4 of the Directive 93/42/EEC on Medical Devices.

We herewith declare that the above mentioned products meet the transposition into national law, the provisions of the following EC Council Directives and Standards. We are exclusively responsible for this Doc.All supporting documentations are retained under the premises of the manufacturer.

DIRECTIVES

General applicable directives:

Medical Device Directive: COUNCIL DIRECTIVE 93/42/EEC concerning medical devices (MDD 93/42/EEC).

Notified Body: TÜV SÜD Product Service GmbH, Ridlerstr. 65, 80339 München, Germany

Identification number: CE0123

(EC) Certificate(s): G1171163196017

Expire date of the Certificate: 2023-04-10

Signature /Date: _____

Name: Mr Lifeng
Position: GM

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50

X

Heal Force

Heal Force International Trading (Shanghai) Co., Ltd6788 Songze Avenue, Qingpu District, Shanghai 201706, China
Tel: 0086-21-62728646 Fax: 0086-21-62710529 Website: www.healforce.com

Authorization Letter

Issuing Date: Dec 10th, 2020

Expiring Date: Dec 10rd, 2021

We, **Heal Force International Trading Shanghai Co., Ltd**, having our office at 6788 Songze Avenue, Qingpu District, Shanghai 201706, China, on behalf of **Shenzhen Hexin Zondan Medical Equipment Co., Ltd**, which is subsidiary companies of **Heal Force International Trading Shanghai Co., Ltd**, for **Heal Force Pulse Oximeter A3** respectively do hereby authorize:

SHIBUMI s.r.l.

Via Carducci 22, 20123 Milano, Italy,

as Heal Force **non-exclusive** distributor of Heal Force products. And **SHIBUMI s.r.l.** has the right to act as Heal Force local agent for registration, sales and services within the authorized territory during the period of validity according to the general business terms and conditions.

Heal Force International Trading (Shanghai) Co., Ltd retains all the right to terminate the authorization by sending **SHIBUMI s.r.l.** written notification 1 month in advance.

Liye Wang

International Sales Director

Heal Force International Trading (Shanghai) Co., Ltd

Address: 6788 Songze Avenue, Qingpu District, Shanghai 201706, China Website: www.healforce.com
Tel: +86 21 6272 8646 Fax: +86 21 6271 0529 E-mail: export@healforce.com

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Product Service

EC Certificate

Full Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)
(Devices in Class IIa, IIb or III)

No. G1 17 11 63196 017

Manufacturer: Shenzhen Hexin ZONDAN Medical Equipment Co., Ltd.

Floor 14, Block D
Dianlian Technology Building
the Crossing between South Circle Road and South Full Road,
Guangming District
518106 Shenzhen
PEOPLE'S REPUBLIC OF CHINA

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)

Eiffestraße 80
20537 Hamburg
GERMANY

Product Category(ies): Pulse Oximeter, Patient Monitor, Fetal Dopplers, Fetal Monitor

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for design, manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex II. This quality assurance system conforms to the requirements of this Directive and is subject to periodical surveillance. For marketing of class III devices an additional Annex II (4) certificate is mandatory. See also notes overleaf.

Report No.: SH1722301

Valid from: 2018-04-11
Valid until: 2023-04-10

Date, 2018-04-11

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TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123

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TÜV SÜD Product Service GmbH · Zertifizierstelle · Ridlerstraße 65 · 80339 München · Germany